LA GESTIONE DELLE MARSI (Medical Adhesive Related Skin Injuries)

MARSI MANAGEMENT (Medical Adhesive Related Skin Injuries)

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KEYWORDS

skin irritations,
erythemas, dialysis,
Marsi, advanced
dressings, adhesives

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ENG

Introduction

Our skin is alive and constantly being renewed. It carries out various highly important functions for the human body, such as protection, temperature regulation, electrolyte balance, metabolism and tactile perception.

Therefore, maintaining skin integrity has always been considered a quality indicator and, when this is absent, patients experience pain and are more vulnerable to infections.

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Skin-related problems experienced by patients subjected to frequent treatments have led scientific research to focus its attention to those referred to as Medical Adhesive-Related Skin Injuries (MARSIs). These are unknown complications, which occur in some patients, requiring frequent medication and repeating dressings using devices with an adhesive component. They are defined as “situations when an erythema and any other type of abnormal skin condition persists for over 30 minutes after removing a plaster or adhesive”.

Substantial clinical evidence proves that, in elderly patients, the incidence of these problems is higher. Furthermore, it has been shown that patients subjected to therapeutic treatments involving changes in normal metabolic parameters are also more susceptible to this type of irritation phenomena.

Research in recent years has classified three types of MARSIs: mechanical, when the corneal layer of the skin is removed after removing the adhesive, erythematous, when non-allergic dermatitis occurs after contact or allergic dermatitis occurs with a cell-mediated response and, finally, other MARSIs, which lead to phenomena such as skin maceration or folliculitis.


In all of these cases, the surface layer of the skin is removed along with the adhesive and this causes damage, not only to the integrity of the skin itself, compromising its barrier function, but also a painful reaction in patients, an increased risk of exposure to infections and delayed healing, resulting in a decreased quality of life for the patient.

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
Mechanical damage with damage to one or more corneal layers, resulting from adhesive removal

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*Dermatitis due to
a cell-mediated
immunological
response after
contact with the
adhesive area*

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*Folliculitis:
inflammatory
reaction caused by
prolonged adhesion
of the adhesive to
the skin*

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The pathophysiology of MARSIs is only partially understood. Damage occurs when resistance to the adhesive exceeds the cohesive strength between skin cells. Consequently, the epidermal layers separate or the epidermis separates completely from the dermis. Even when there is no visible trauma, common adhesive removal causes the detachment of the epidermal surface layers and repeated application and removal compromises the skin barrier function, resulting in inflammation and a slowed healing response of the wound.

In more severe cases, adhesives can even cause damage to tissues underlying the skin.

Although this type of condition appears to be mild, MARSIs management, in terms of cost and use of nursing staff, can be extremely costly.

These reactions often occur despite the use of correct dressing removal technique by medical staff, leading to suggest that nothing can be done to resolve this problem. In fact, it is necessary to emphasise that this happens when the bond between the adhesive and the surface layer of the skin is stronger than that created between the corneal layer and the underlying layers.

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The composition of the adhesive itself is therefore crucial for preventing these problems. Traditional adhesives, especially acrylic ones, increase the strength of adhesion to the skin, as they are obtained from aqueous dispersion acrylic resins. During the coating stage, water is evaporated. They are characterised by good resistance to heat, medium adhesion, good resistance to solvents and ageing, but also poor resistance to humidity. To date, acrylic materials have, however, been identified

as an important source of contact sensitisation, and are responsible for many allergy phenomena due to their specific chemical composition. With this evidence, Emodial research has focused on new solutions.

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The Emodial Experience

In patients undergoing dialysis therapy, which requires continuous and prolonged dressings, MARSI management has become an urgent requirement.

In patients with chronic kidney disease, substances that should normally be eliminated from the body accumulate, causing progressive metabolic abnormalities in many organs and systems, leading to many dermatological complications, especially in patients undergoing dialysis. The most common dermatological complication is skin xerosis, known as uremic xerosis, a condition that aggravates persistent "uremic pruritus", especially in the legs, trunk and hands.

Uremic xerosis is a chronic condition characterised by the presence, on the body, of dry, rough and

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squamous skin, the cause of which is due to the deficiency of the hydro-electrolytic balance that is typical in patients with kidney failure and skin glycerol deficiency, resulting from the reduced activity of the sebaceous glands.

Reduced hydration of the skin and deficiency of the main water transport molecule, glycerol, result in intense skin dehydration.


These types of patients are particularly exposed to skin disorders. This evidence has made it necessary to seek an alternative to alleviate or resolve these skin reactions and thus improve quality of life that has been compromised by this condition.

A viable solution for preventing MARSIs was identified in the use of hypoallergenic adhesives that do not damage the skin and therefore do not have the features typical of normal medicinal adhesives which, in the long term, cause irritation.

SoftTouch dressings are coated in a layer of hydrophobic silicone, a synthetic material comprising a long polymer chain that consists of shorter, repeated silicone chains, combined with carbon, hydrogen and oxygen.

The chemical properties of this substance make it particularly adaptable and conformable and enable it to adhere well to dry surfaces. Silicones have low toxicity, are not associated with adverse reactions and cannot be absorbed by the body. This makes them ideal for use in dressings.

SoftTouch dressings do not stick to the wound bed, but adhere to the skin surrounding it and are specifically designed to minimise trauma due to

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removal and do not leave residues on the skin.


They are "tack" dressings, that is, dressings that can immediately adhere to the surface and, unlike acrylates, do not increase in firmness over time. This feature of silicones is related to their ability to create many points of contact between the dressing itself and the surface layer of the skin, generating the formation of a close connection with the wound, which prevents skin maceration, thus preventing the release of exudate.

They are also ideal dressings as they enable adhesion to be maintained throughout the duration of their action and enable removal without trauma to the wound bed.

*Emodial **SoftTouch** Line dressings comprise a delicate silicone adhesive that is supported by a layer of unwoven fabric and specifically indicated for wide range of exudative lesions, including ulcers, pressure ulcers, diabetic foot, second-degree burns, post-operative lesions and skin abrasions.*

When in contact with the exudate, the silicone adhesive layer conforms to the wound bed to ensure optimal absorption, by minimising dead spaces and absorbs the exudate vertically, with a high retention capacity, thus reducing the risk of damage to the damaged skin caused by maceration. It also provides secure adhesion that enables the dressing to be kept in position, ensuring, at the same time, an almost painless removal.

The soft and flexible structure of the dressing ensures perfect conformability to the lesion and to the body, making the dressing comfortable to use. This adhesive allows the dressing to adhere to the skin again if it is repositioned.

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
Clinical Evidence

*Emodial, in partnership with some of the most important hospital centres specialising in dialysis, recently conducted a clinical survey aimed at testing **SoftTouch** adhesive technology on a sample of patients normally treated with traditional adhesives.*

The survey involved patients filling in a questionnaire comprising a series of questions to assess the effects of these dressings.

They were requested to give their opinion in terms of:

- Adhesiveness
- Trauma caused by removal
- Skin irritation
- Comfort

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The results were very encouraging: patients declared a significant improvement in aspects taken into consideration.

*Emodial **SoftTouch** line dressings were tested on 30 patients from various Italian hospitals. The graph shows the results of the assessment on devices with silicone adhesive, compared with a previously used standard dressing.*

Bibliografia

Bibliography

Acton C (2007) The holistic management of chronic wound pain. *Wounds UK*. 3, 1, 61-69.

Cassidy C et al (2005) Biobane versus duoderm for the treatment of intermediate thickness burns in children: a prospective, randomised trial. *Burns*. 31, 7, 890-893. European Wound Management Association (2002)

Gotschall CS et al (1998) Prospective, randomised study of the efficacy of Mepitel on children with partialthickness scalds. *Journal of Burn Care Rehabilitation*. 19, 279-283.

Harms-Ringdahl K et al (1986) Pain assessment with different intensity scales in response to loading of joint structures. *Pain*. 27, 3, 401-411.

Hollinworth H, Collier M (2000) Nurses' views about pain and trauma at dressing changes: results of a national survey. *Journal of Wound Care*. 9, 8, 369-373.

Keck JF et al (1996) Reliability and validity of the Faces and Word Descriptor Scales to measure procedural pain. *Journal of Pediatric Nursing*. 11, 6, 368-374.

Platt AJ et al (1996) A comparative study of silicone net dressing and paraffin gauze dressing in skin-grafted sites. *Burns*. 22, 7, 543-545.

Soon K, Acton C (2006) Pain-induced stress: a barrier to wound healing. *Wounds UK*. 2, 4, 92-101.

Thomas S (2003) Atraumatic dressings. *World Wide Wounds*. Available online at: www.worldwidewounds.com.

White R (2005) Evidence for atraumatic soft silicone wound dressing use. *Wounds UK*. 1, 3, 104-109.