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ALGINATE BASED DRESSINGS FOR ARTERIOVENOUS FISTULA MANAGEMENT

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KEYWORDS

haemostasis, advanced dressings, dialysis, fistula, coagulation, exudates, healing, coagulopathies

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ENG

Emodial's experience alongside the main dialysis centres in Italy and abroad has revealed an increased presence of patients with haemostasis problems or with specific coagulation-related diseases. This condition, along with nephrological disease, has driven the search for dressings capable of limiting the consequences of the two diseases which, acting simultaneously on the patient, often lead to a particularly delicate clinical picture.

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Patients with haemostasis problems and coagulopathies comprise a large part of the population. In Italy, for example, referring only to patients undergoing anticoagulant therapy, a prevalence of 911,000 cases is estimated, with a percentage of 1.58% of the total population. Many elderly patients, diabetics, certain patients with vitamin K deficiency and some patients treated with oral anticoagulants (such as Coumadin) or heparin can be included in this category.

Haemodialysis itself uses heparin to protect the bloodstream from the risk of clotting. Abnormalities in the lumen and intravessel wall of the arteriovenous fistula, respectively, the presence of stenosis and aneurysms, create haemostasis problems, regardless of the presence of coagulopathies or controlled heparin administration.

Serious liver diseases, such as cirrhosis and fulminant hepatitis, can also interfere with correct haemostasis by compromising the synthesis of coagulation factors, normally synthesised in the liver.

A large number of patients suffer from hereditary forms. The most common disorder is haemophilia, a typically male disease caused by a factor VIII coagulation deficiency (haemophilia A) or factor IX coagulation deficiency (haemophilia B). Other hereditary coagulopathies are Von Willebrand's disease and some forms of platelet disease due to intrinsic platelet function abnormalities.

For this category of patients, research has identified advanced calcium alginate-based dressings, a viable ally for resolving haemostasis problems. Due to its properties, calcium alginate is a substance that can be used in a preventive manner, extending its benefits associated with reduced bleeding times to patients not suffering from haemostasis problems.

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
What is alginate?

Alginates are alginic acid, or, more simply, algin salts. Algin is a polymer derived from the cell wall of some brown algae, such as laminaria and fucus. They are a family of polysaccharides, characterised by a very variable chemical composition, molecular weight and functional properties. Chemically, they are linear polymers of β -D-mannuronic (M) acid and α -L-guluronic (G) acid, which exhibit a different structural geometry. This difference is due to the ability of bivalent cations, such as Ca^{2+} , to bind to alginate molecules, only at the G blocks, causing them to form a gel.

Alginate is soluble in aqueous solution and capable of changing its viscosity depending on the interactions in progress. It is therefore capable of creating stable three-dimensional structures in an aqueous environment and creating intermolecular bonds by chelation of bivalent cations.

This feature of forming a gel by interacting with an aqueous environment, on the contact surface, has always been exploited in the field of wound care.

Alginate interacts with the exudate of the lesion and forms a soft gel which keeps the healing environment moist. Comprising cation ions, it may promote the coagulation process. The high absorbent potential of this category of dressings and their ability to adapt to

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any base, facilitates the removal of bacterial residues through gelling in the gelified matrix.

Alginates can be used as primary dressings in the case of draining lesions of a partial or total thickness, moderate lesions with abundant exudate, tunnel lesions and infected lesions.

Those containing calcium ions can be indicated in bleeding lesions, as they promote haemostasis. When in contact with blood, gelified calcium alginate, causing calcium ions to pass into the wound, can promote one of the first stages in haemostasis, platelet aggregation and can shorten the time needed for coagulation.

By assuming a gelatinous consistency when in contact with blood or with exudates of a weeping lesion, the calcium alginate creates a sort of artificial skin that protects the wound from external impacts. This gel also maintains a certain rate of humidity in the wound, accelerating the tissue repair process.

Generally, alginate-based dressings are extremely adaptable to any type of surface on lesions of various sizes, shapes and depths. They conform to the bed of lesion and thus promote healing.

They are highly absorbent (they can absorb exudate equal to 10-20 times their weight) and can minimise loss of exudate that can damage the surrounding skin. They are very easy to remove as the gel formation prevents the dressing from sticking to newly-formed tissue.

They do not leave residues in the wound and this promotes a more comfortable dressing change for patients. They can also be removed in a single piece.

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Alginate in Dialysis: Emodial Dressing

Emodial alginate-based medical devices comprise absorbent pads made from calcium alginate fibre on various types of support, such as a cotton bandage. In the presence of exudates or other body fluids containing sodium ions, the fibres absorb liquids and the calcium ions present are partially replaced by sodium ions, leading to the formation of a gel. This provides a micro-environment that facilitates haemostasis and healing of the wound.

These devices are intended for wound management, including chronic wounds such as ulcers and acute wounds such as abrasions, lacerations and post-surgical wounds. They are also intended for managing arteriovenous fistulas (AVF) with haemostasis problems.

*The structure of the multi-layer pad on the **BETAFIX[®] ALGINATE AG** device, an AV fistula bandage with a calcium alginate-based pad, makes it particularly suitable for use in severe bleeding, thanks to the synergistic action performed by the alginate and the mechanical action, performed by the overall action of the bandage.*

Calcium ions associated with the dressing also when they are released into the wound through the blood vessels trigger the formation of a clot and accelerate haemostasis and activate the key cells in the scarring process, known as macrophages. Also, when applied directly onto the wound, these medical devices compress the wound and act as a mechanical barrier to protect it from the external environment, creating a humid environment and, therefore, promoting tissue regeneration.

85.6% of patients undergoing haemodialysis with coagulation problems have shown a reduced bleeding time of 5 minutes.

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Emodial conducted a pre-clinical study on voluntary patients undergoing extracorporeal dialysis with an arteriovenous fistula, in various Italian hospitals, by asking doctors to fill in a questionnaire on the efficacy of the calcium alginate-based bandage.

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Doctors were asked to compare this new device with that normally used, in terms of:

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- *amount of blood released*
- *speed of haemostasis*
- *skin irritation*
- *scarring*

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The results gathered were very encouraging, as can be seen in the table below.

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Outcome of the reduced haemostasis time assessment with a calcium alginate-based pad vs. standard dressing

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
The final line shows that 11 out of 13 patients who had a coagulation time greater than 21 minutes has a reduction in bleeding time to less than 5 minutes and 2/13 less than 10 minutes: no patient took more than 10 minutes to reach haemostasis.

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Reduced AVF haemostasis times following application of a calcium alginate-based dressing

REFERENCE RECORDS

- 1 THOMAS, S.
Alginate dressings in surgery and wound management - Part 1.
Journal of wound care, 2000, 9.2: 56-60.
- 2 THOMAS, S.
Alginate dressings in surgery and wound management - Part 3.
Journal of wound care, 2000, 9.4: 163-166.
- 3 CHABY, Guillaume, et al.
Dressings for acute and chronic wounds: a systematic review.
Archives of dermatology, 2007, 143.10: 1297-1304.
- 4 QIN, Yimin.
Alginate fibres: an overview of the production processes and applications in wound management.
Polymer international, 2008, 57.2: 171-180.
- 5 JUDE, E. B., et al.
Prospective randomized controlled study of Hydrofiber® dressing containing ionic silver or calcium alginate dressings in non-ischaemic diabetic foot ulcers.
Diabetic Medicine, 2007, 24.3: 280-288.
- 6 A. GREN, M. S.
Four alginate dressings in the treatment of partial thickness wounds: a comparative experimental study.
British journal of plastic surgery, 1996, 49.2: 129-134.
- 7 BELMIN, Joël, et al.
Sequential treatment with calcium alginate dressings and hydrocolloid dressings accelerates pressure ulcer healing in older subjects: a multicenter randomized trial of sequential versus nonsequential treatment with hydrocolloid dressings alone.
Journal of the American Geriatrics Society, 2002, 50.2: 269-274.
- 8 DOYLE, James W., et al.
Effect of calcium alginate on cellular wound healing processes modeled in vitro.
Journal of biomedical materials research, 1996, 32.4: 561-568.
- 9 KNEAFSEY, B.; O'SHAUGHNESSY, M.; CONDON, K. C.
The use of calcium alginate dressings in deep hand burns.
Burns, 1996, 22.1: 40-43.